

Ricordi Chamber O-Ring Model No. OR-XXX



Each Ricordi® Chamber uses one silicone o-ring. They can be sterilized with steam autoclave or EtO. To ensure proper sealing of the chamber, make sure there are no cracks or cuts on the o-ring.

- Temperature Range : -65° to +450°F
- Durometer: Shore A 70
- FDA Specification: FDA Compliant
- Color: red-orange

To order, please specify the chamber model.

Hole size μm	Percentage of open area
50mL	OR-50
125mL	OR-125
250mL	OR-250
400mL	OR-400
500mL	OR-500
600mL	OR-600
700mL	OR-700
1000mL	OR-1000